DELABIE

BIOFIL

NT30250

Indice 1

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1 REFERENCES

- BIOFIL Shower head water filters: 30261.
- BIOFIL Tap water filters with straight flow outlet: 30251.
- BIOFIL Tap water filters with shower flow outlet: 30250.

To guarantee the optimal performance of the BIOFIL disposable water filters, please ensure that you follow the instructions for installation (see 5) and read the warnings carefully (see 6).

2 GENERAL SPECIFICATIONS

- Disposable water filters.
- Effective against bacteria including opportunistic waterborne pathogenic microorganisms.
- 2 months.

BIOFIL filters are **class I medical devices** in Europe.

3 INDICATIONS FOR USE

BIOFIL 2-Month water filters provide a barrier against waterborne microorganisms such as bacteria (*Legionella pneumophila*, *Pseudomonas aeruginosa* and other Gram negative bacilli, nontuberculous Mycobacteria), fungi (*Aspergillus sp*, *Fusarium sp*), protozoa and particles originating in the water distribution system. The microfiltration membrane protects **immune-suppressed patients** and vulnerable people. Filters are necessary for the care of hospitalised, **immune-suppressed patients** in high-risk areas such as bone marrow transplant and organ transplant units, oncology services, intensive care units, resuscitation, neonatal services and any other places that care for this type of patient.

4 DESCRIPTION

The BIOFIL filters are intended to produce bacteriologically controlled water at the point-of-use. These devices retain bacteria, including opportunistic waterborne pathogenic microorganisms such as *Legionella spp.*, *Legionella pneumophila*, *Pseudomonas aeruginosa* and others via a microfiltration membrane. This microfiltration technology ensures a higher bacteriological quality than the potable water of the water plumbing system. Please note, the BIOFIL devices do not contain latex, nor tissue or by-products of animal origin, neither are they derived from blood or any substance considered as a medicine.

BIOFIL disposable water filters are designed to ensure water safety at the point-of-use. They are recommended in the procedures for **management and prevention of healthcare associated infections**.

BIOFIL filters use a 0.1µm nominal / 0.2µm absolute tubular membrane filtration technology of sterilising grade (see 9-A). The filters provide a **high-performance level of water microfiltration** while delivering a **comfortable flow** for users.

Microbiological challenges, conducted in laboratories with *Brevundimonas diminuta*, *Legionella pneumophila* and *Pseudomonas aeruginosa* species on BIOFIL filters, demonstrate a **retention efficiency of 7 log** (see 9-2) for all water-borne microorganisms greater than 0.1µm nominal / 0.2µm absolute and qualify the filter as a sterilising grade filter (see 9-A).

5 INSTALLATION INSTRUCTIONS

Warnings:

- Filter installation must conform with professional aseptic techniques. Hands must be disinfected prior to installation. Inappropriate installation may result in a source of contamination or retro-contamination.
- Make sure that the individual unit packaging (plastic bag) is not punctured or torn and that the filter is not damaged. If this is the case, do not use the device, as the performance can no longer be ensured.

General information:

- The threaded shower head filter: can be attached quickly, screw in clockwise until the seal is right up against the flexible hose. The supplied seal must be fitted to ensure waterproofing and prevent possible leaks. It is advisable to replace the seal each time the shower head is changed.
- The clip-on tap filters: must be fitted using the correct push-fit connectors enabling the filter to be connected and disconnected with ease to the water control. Contact us for the list of available connectors.

Installation:

1. Turn off the water supply.
2. If this is the first installation: remove the original shower head or fit the adaptor, if necessary, for the tap, wall-mounted shower (see removal below).
3. Remove the used BIOFIL filter.
4. Carefully open the individual unit packaging (bag), making sure that the filter remains inside the bag. **Do not touch the end of the tap or the shower head with your fingers to avoid contamination.**
- 5 **Keeping the filter in the bag**, screw the shower head filter onto the end of the hose or clip on the tap or wall-mounted shower filter.
6. The bag can then be removed.
7. Turn the water on for a few seconds to check that the adaptor(s) and the filter are **perfectly watertight**. In the case of leaks, repeat the installation procedure, using a new filter if needed.

Removal:

- 1. Turn off the water supply.
- 2. Remove the used shower filter by unscrewing it from the end of the flexible hose or unclip the shower filter or tap filter from the adapter connector.
- 3. Dispose of the used filter according to your institution's established procedures.

Use in a non-controlled environment can be a source of retro-contamination. Any filter that is believed to be contaminated must be replaced, especially if the contaminant is visible.

Waste management:

Used filters do not present any special infectious risk. They are classified as **normal household waste**. The installation guide and the white carton can be recycled (paper/card).

Maintenance:

To facilitate the tracking of maintenance procedures, a removable adhesive label (located on the filtration device) can be used to record important information for your safety records (date of manufacture, date of removal, batch number and filter unique identification number, etc.). Stick on to the tracking page and write down the installation and removal date. This information is also available as the UDI (Unique Device Identification) which is a data matrix with information in the form of (01) GTIN: Global Trade Item Number, (10) Lot number, (17) expiration date, (21) Unit Number, (240) DELABIE part number.

Manufacturing date, Expiration date and Lot number are also clearly identified on device labels and packaging. The lifespan is identified by a numerical symbol featuring the number of months.

6 CONTRAINDICATIONS

BIOFIL Disposable Water Filters:

- Must not, under any circumstances, be used to produce water for **pharmaceutical purposes** (i.e., water that will be used in solutions intended to be administered orally, injected, or perfused).
- Are single use (disposable) or re-installed on another point-of-use.
- Have a **limited lifespan** according to the filter used, and an **expiration date** indicated on the label. These must be adhered to.
- Must not come into direct contact with germ sources (dirty skin or hands, contaminated aerosols, body fluids, etc.).
- Must not be unscrewed or disconnected during the life of the filter except if the water distribution system is disinfected. In this instance, replace with a new filter.
- Must not be struck or dropped (fall). Filter performance is not guaranteed if the filter casing is physically damaged. Do not attempt to open the device.
- Shower head filters must always be replaced in the wall-mounted bracket after use.
- All filters must be regularly cleaned then disinfected according to your institution's established procedures, starting at the water outlet, using a clean disinfectant wipe.

7 TECHNICAL SPECIFICATIONS

Conditions for use:

The performance of BIOFIL disposable water filters is guaranteed when the water pressure and temperature remain between 1 to 5 bar and 5°C and 42°C, respectively. Flow rates measured under normal conditions are as follows:

BIOFIL filter					
	Water pressure (bar) static or dynamic				
	1	2	3	4	5
Flow rate (l/min) - Tap 2M	3	4.4	5.5	6.5	7.4
Flow rate (l/min) - Shower Head 2M	6	8.8	11.2	12.8	13

Heat and chemical resistance:

BIOFIL disposable water filters are resistant to the thermal and chemical shocks that frequently occur during disinfection of the water supply system (see 6-1). They remain effective after treatment at **70°C for 30 min** (cumulative maximum) or a dose of active chlorine at **100 ppm for 60 min** (cumulative monthly maximum): BIOFIL 2-month filters: 2 hours at 100ppm.

Blockage tolerance:

The flow rate at the point of use may decrease somewhat before the end expiration date if there are unusually high levels of particles or suspended microorganisms in the water supply (for example, after a maintenance procedure has been performed). If the flow rate becomes too slow for user comfort, simply replace the filter.

8 PRODUCT COMPLIANCE

BIOFIL Disposable Water Filters:

- Are **Class I Medical Devices in France** - complying with European directive 2017/745. **CE marking: year 2021.**
- Comply with drinking water certifications (**ACS-France, KTW-Germany, WRAS-United Kingdom**).

BIOFIL Disposable Water Filters are manufactured in France by:
DELABIE | 18, rue du Maréchal Foch | 80130 FRIVILLE | France
For more information, visit our website at: www.delabie.co.uk
or contact us: **tel.:** +33 (0)3 22 60 22 74 - **email:** sav@delabie.fr

References:

- A. Sterilising grade filter is defined in “American standard Test Method ASTM F838-20 and in FDA - Guidance for Industry - Sterile Drug Products Produced by Aseptic Processing - Current Good Manufacturing Practice” as a filter ability to sustain a bacterial challenge of 10⁷ CFU of *Brevundimonas diminuta* per cm² of filtration surface.

Product label - meaning of symbols

		
Catalogue number	Batch code	Use by
		
Sterilised using irradiation	Do not resterilise	Do not use if packaging is damaged
		
Free from natural latex	Manufacturer	Compliant with the general requirements for safety and performance
		
Non-sterile	Consult instructions for use	

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